

■ C120030

PATIENT STRETCHER

- Epoxy painted steel frame, platform fitted with enamelled steel weld mesh. Backrest easily adjustable.
 - Supplied with I.V. stand and side rails.
 - Mounted on four \varnothing 200 mm casters, two of which with brakes.
 - Dimensions: 75 x 198,5 x 75 cm. Capacity: 200Kg
-

BARELLA

- Struttura in metallo, rivestimento epossidico.
- Base rigida in rete metallica. Alzatesta.
- Fornito con asta portaflebo e spondine laterali
- Quattro ruote \varnothing 200 mm, di cui due con freno indipendente.
- Dimensioni: 75 x 198,5 x 75 cm. Capacità: 200Kg