

## PATIENT STRETCHER

- **C120031** Patient stretcher with variable height
- **C120032** Patient stretcher variable height with Trendelenburg.
- Main frame made of round section steel tube, platform fitted with removable enamelled steel weld mesh.
- Supplied with I.V. stand and side rails.
- Backrest easily adjustable by gas pump. I.V. stand.
- Mounted on four ø 200 mm. casters, two with brakes.
- Dimensions: 65 x 197 x 50cm
- Variable height 50/89cm. Capacity: 200Kg

## BARELLA

- **C120031** Fornita con altezza variabile
- **C120032** Fornita con altezza variabile e posizione Trendelenburg (+7°)
- Robusta struttura in metallo, rivestimento epossidico.
- Base rigida in rete metallica, 2 sezioni..
- Spondine laterali, asta portaflebo, supporto per bombola di ossigeno
- Alzatesta regolabile con pompa a gas.
- Quattro ruote ø 200 mm, di cui due con freno indipendente.
- Dimensioni: 65 x 197 x 50cm
- Altezza variabile 50/89cm. Capacità: 200Kg