

■ C120040

EMERGENCY STRETCHER

- Oleo-dynamic semi-automatic pedal-controlled stretcher
- Steel frame, bed board and side rails by PP plastic
- Centralized brake ensures perfect stability and hold
- Variable height and TR/RTR positions (0-13°) adjustable by pedal
- Supplied with: two folding siderails, adjustable I.V. pole, an oxygen-tank holder and a wide tray
- Mounted on four swiveling 150 mm castors and 5th central directional wheel
- The stretcher is supplied with flame retardant mattress, density 30 kg/m³ with washable and waterproof cover
- Dimensions: 61 x 190,5 x 58/81,5 cm. Capacity: 200 Kg

BACKREST
ALZATESTA

VARIABLE HEIGHT
ALTEZZA VARIABLE

TRENDELENBURG
and REVERSE

BARELLA D'EMERGENZA

- Barella oleo-dinamica semi-automatica con comando a pedale
- Struttura in acciaio, piano di coricamento e sponde laterali in plastica PP
- Il freno simultaneo garantisce stabilità e tenuta perfette
- Altezza variabile e posizioni TR / RTR (0-13 °) regolabili tramite pedale
- Fornito con: due sponde pieghevoli, Asta porta flebo regolabile, supporto per bombola d'ossigeno e un ampio vassoio
- Montato su quattro ruote orientabili da 150 mm e una quinta ruota direzionale
- La barella è fornita di materasso ignifugo, densità 30 kg / m³ con fodera lavabile e impermeabile
- Dimensioni: 61 x 190,5 x 58 / 81,5 cm. Portata: 200 Kg