

■ C120039

EMERGENCY STRETCHER

- Manual transfer stretcher single crank with height adjustment
- Main frame made of powder-coated metal sheet, bed board and side rails made of PP materials
- Mounted on \varnothing 150 mm castors with central locking system and 5th direction wheel under the chassis
- Size of bed board: 190 x 61 cm, back-rest lifting angle: 0-60° +/- 5°
- Drop-down side rails, standard accessories: mattress, I.V. pole
- The stretcher is supplied with flame retardant mattress, density 30 kg/m³ with washable and waterproof cover
- Dimensions: 202 x 74 x 58/88 h cm. Capacity: 200 Kg

BARELLA D'EMERGENZA

- Barella manuale con altezza variabile regolabile tramite manovella
- Struttura in acciaio, piano di coricamento e sponde laterali in plastica PP
- Montato su quattro ruote orientabili da 150 mm e una quinta ruota direzionale
- Fornito con: due sponde laterali pieghevoli, Asta porta flebo regolabile, supporto per bombola d'ossigeno
- Dimensioni del piano: 190 x 61 cm, schienale: 0-60° +/- 5°
- La barella è fornita con materasso ignifugo, densità 30 kg/m³ con fodera lavabile e impermeabile
- Dimensioni: 202 x 74 x 58/88 h cm. Portata: 200 Kg