

■ 01.6265

PROCEDURE CART

- Made of aluminum-alloy, ABS plastic and stainless steel.
- 6 layers of drawers which can be pulled out from both sides, self-locked when closed.
- 50 compartments card slot holder at the front of every compartment
- 2 waste bins, 1 transparent file box, 1 sliding table and 1 hanging device supplied with 5 hooks.
- 4 castors diam. 10 cm, 2 with brakes
- Dimensions: 85 x 53 x 101 cm

CARRELLO PER MEDICAZIONI

- Realizzato in lega di alluminio, plastica ABS e acciaio inossidabile.
- 6 cassetti estraibili da entrambi i lati, autobloccanti quando chiusi.
- 50 scomparti di supporto per slot per schede nella parte anteriore di ogni scomparto
- 2 pattumiere, 1 scatola per limatura trasparente, 1 tavolo scorrevole e 1 appendino con 5 ganci.
- 4 ruote diametro 10 cm, 2 con freno.
- Dimensions: 85 x 53 x 101 cm